



# ORGANICS

---

Arizona Department of Health Services

Mary Cornell

# Method Updates

- ADHS will no longer license method 8141A. Method 8141B will be offered instead. <sup>14</sup>
- Both 8260B and 8260C will be offered as well as both 8270C and 8270D.
- The version of method 8000 you will reference will depend on which determinative method you will be using. <sup>15</sup>

Method 8000	Determinative Methods
8000C	8141B 8260B 8270C
8000D	8260C <sup>16</sup> 8270D

## Slide 2

---

- I4 May director approve 8141A if permit requires it...check with Steve.  
labuser, 9/2/2016
- I5 May not pick and choose requirements/allowances from 8000C and 8000D. Must follow one or the other.  
labuser, 9/2/2016
- I6 8000D not explicitly assigned to 8141B in the draft rules. This was our mistake. Too late to correct. We will be addressing this.  
labuser, 9/2/2016

# General SW-846 Method Updates

- Improved overall method formatting for consistency with SW-846 methods style guidance.
- Formatting updated to Microsoft Word.docx
- Many minor editorial and technical revisions made throughout to improve method clarity. 12

### Slide 3

---

I2

This presentation will not cover these minor revisions in detail. Only major revisions will be covered.

labuser, 9/9/2016



## Method 8000D

# Determinative Chromatographic Separations

Revision 4, July 2014

# 8000D Updates

**Carrier Gas** – Optional use of hydrogen gas for GC to address helium shortages (Section 1.3.1)

## **Initial Demonstration of Proficiency (IDP)**

- Replaces “Initial Demonstration of Capability (IDC)” in 8000C
- Section 9.3.3 – Additional text to indicate that samples used for IDP analyses should be from the same source as the calibration standard.
  - Now consistent with the language in methods 8260 and 8270 (Appendix A)
- Separate language for preparation/extraction chemists and instrument chemists.
  - This allows the chemist to only perform an IDP for the procedure that they actually do. (Appendix A)

## 8000D (cont.)

- **MQL vs. LLOQ** - Method 8000D will now use the terminology Lower Limit of Quantitation (LLOQ) in place of Method Quantitation Limit (MQL).
- **Relative Standard Error (RSE)** - Included in Section 11.5.4.2 as an alternative for the determination of calibration acceptability.
  - RSE acceptance criteria is the same as the RSD limit for RF in the determinative method.
  - If the RSD limit is not defined in the determinative method, the limit should be set at  $\leq 20\%$  for good performing compounds and  $\leq 30\%$  for poor performing compounds.



# 8000D Calibration Verification

- LCS may serve as the CCV (Section 9.4.1)
- For multi-analyte methods, you will not have to reanalyze if  $\leq 10\%$  of compounds fall outside of the acceptance criteria. (Section 11.7.1)
  - Any detected analytes exceeding the limit **must be reported as estimated**.
  - If the analyte recovers below the acceptance criteria and the analyte is non-detect in the sample, a sensitivity verification standard  $\leq$  LLOQ should be analyzed. The analyte should be detected in the LLOQ standard and meet qualitative identification criteria (e.g. RT, qualifier ions, signal-to-noise, etc.).



## Method 8260C

# Volatile Organic Compounds by Gas Chromatography/Mass Spectrometry

Revision 3, August 2006

# 8260 ICAL

8260B	8260C
Calibration Check Compounds (CCC)	No CCCs
System Performance Check Compounds (SPCC)	No SPCCs
%RSD $\leq 15$ for all target compounds; %RSD for CCCs $\leq 30$ (Section 7.3.6.3)	%RSD criteria for all targets $\leq 20$ (Section 11.3.4.1)
No minimum RF requirement for target analytes	Recommended minimum RF values specified in Table 4 (Section 11.3.4.1)
No re-fitting of calibration standards required	If linear calibration used must verify the RL by re-calculating concentrations in the lowest calibration standard using calibration curve. Acceptance criteria is 70-130% Rec. (Section 11.4.5.6)
Must recalibrate if 1 or more CCCs out or $>20\%$ of targets exceed 30% RSD (Section 7.3.6.4)	Must recalibrate if $>10\%$ of targets exceed the %RSD or regression criteria (Section 11.3.4.2)

## 8260C ICAL – Additional Requirements

18

- May use a single surrogate concentration in ICAL if autosampler is used (Section 11.3.3 Note)
- SIM - must monitor a minimum of 2 ions per analyte, primary and confirmation (Section 11.5.12)

## Slide 10

---

I8

Most labs already complying with the following. These are just specifically addressed now.

labuser, 9/2/2016

# 8260 Calibration Verification

8260B	8260C
No ICV requirement	Must be a second source standard. Suggested acceptance criteria 70-130% (Section 11.4.2)
No RF requirement for CCVs	Evaluate minimum RFs specified in Table 4 of 8260C (Section 11.4.5.1)
No IS requirement for CCVs	The IS area counts must be -50-200% of area counts in the associated mid-level ICAL std. (Section 11.4.7)
For a CCV, CCC %Drift criteria $\leq 20$ (Section 7.4.5.2)	%D $\leq 20$ for all targets. Recalibrate if $>20\%$ targets exceed the %D criteria. (Section 11.4.5.4) <ul style="list-style-type: none"><li>Failed compounds may be reported as non-detects if there is adequate sensitivity to detect the compound at the quantitation limit.</li><li>Failed compounds that are present in the sample must be reported as estimated values.</li></ul>

# 8260C Additional Updates

- **Section 7.9** – Dibromofluoromethane no longer listed as a recommended surrogate, but may be used depending upon the analysis requirements.
- **Section 7.13** - LCS and MS should be prepared from the same source as the calibration standards.<sup>I1</sup>
- **Section 7.14** - Standards may be stored <6°C (previously<sup>I9</sup> < -10°C) unless otherwise specified by the manufacturer.
- **Section 11.6.1.4** - Sufficient GC resolution is achieved if the height of the valley between two isomer peaks is <50% (previously <25%) of the average of the two peak heights. Otherwise, structural isomers are identified as isomeric pairs.

## Slide 12

---

- I1 Section 7.13 already allowed by Director Approved Method Modification #20.  
labuser, 9/2/2016
- I9 Recommend to continue storing volatile standards <-10 degrees C.  
labuser, 9/2/2016





## Method 8270D

# Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry

Revision 5, July 2014

# 8270 DFTPP Tune Criteria

## 8270C

51	30-60% of mass 198
68	< 2% of mass 69
70	< 2% of mass 69
127	40-60% of mass 198
197	< 1% of mass 198
198	Base peak, 100% relative abundance
199	5-9% of mass 198
275	10-30% of mass 198
365	> 1% of mass 198
441	Present but less than mass 443
442	> 40% of mass 198
443	17-23% of mass 442

## 8270D

<b>10-80%</b> of mass 198
< 2% of mass 69
< 2% of mass 69
<b>10-80%</b> of mass 198
<b>&lt; 2%</b> of mass 198
Base peak, <b>or &gt; 50% of mass 442</b>
5-9% of mass 198
<b>10-60%</b> of mass 198
> 1% of mass 198
Present but <b>&lt; 24% mass 442</b>
<b>Base peak, or &gt; 50% of</b> mass 198
<b>15-24%</b> of mass 442

# 8270 ICAL

8270C	8270D
Calibration Check Compounds (CCC)	No CCCs
System Performance Check Compounds (SPCC)	No SPCCs
%RSD $\leq 15$ for all target compounds; %RSD for CCCs $\leq 30$ (Section 7.3.5.2)	%RSD criteria for all targets $\leq 20$ (Section 11.3.4.1)
No minimum RF for target compounds	Recommended minimum RF values specified in Table 4 (Section 11.3.4.1)
No re-fitting of calibration standards required	If linear calibration used must verify the RL by re-calculating concentrations in the lowest calibration standard using calibration curve. Acceptance criteria is 70-130% Rec. (Section 11.4.5.6)
Must recalibrate if 1 or more CCCs out or >20% of targets exceed 30% RSD (Section 7.3.5.3)	Must recalibrate if >10% of targets exceed the %RSD or regression criteria (Section 11.3.4.2)

Identical to 8260C ICAL updates

# 8270 Calibration Verification

8270C	8270D
No ICV requirement	Must be a second source standard. Suggested acceptance criteria 70-130% (Section 11.4.2)
No RF requirement for CCVs	Evaluate minimum RFs specified in Table 4 of 8260C (Section 11.4.5.1)
No IS requirement for CCVs	The IS area counts must be 50-200% of area counts in the associated mid-level ICAL std. (Section 11.4.7)
For a CCV, CCC %Drift criteria $\leq 20$ (Section 7.4.5.2)  Also identical to 8260C	%D $\leq 20$ for all targets. Recalibrate if $>20\%$ targets exceed the %D criteria. (Section 11.4.5.4) <ul style="list-style-type: none"><li>Failed compounds may be reported as non-detects if there is adequate sensitivity to detect the compound at the quantitation limit.</li><li>Failed compounds that are present in the sample must be reported as estimated values.</li></ul>

## 8270D Additional Updates

- **Section 8.2** - Standards may be stored  $<6^{\circ}\text{C}$  (previously  $< -10^{\circ}\text{C}$ ) unless otherwise specified by the manufacturer.
- **Section 11.6.1.4** - Sufficient GC resolution is achieved if the height of the valley between two isomer peaks is  $<50\%$  (previously  $<25\%$ ) of the average of the two peak heights. Otherwise, structural isomers are identified as isomeric pairs.

## Slide 17

---

**I10** Section 7.13 already allowed by Director Approved Method Modification #20.

labuser, 9/2/2016

**I11** Recommend to continue storing volatile standards <-10 degrees C.

labuser, 9/2/2016



## Method 8141B

# Organophosphorus Compounds by Gas Chromatography

Revision 2, February 2007

# 8141B Updates

- Solid samples may now be extracted by EPA Method 3550 - Ultrasonic Extraction (Section 2.1.2)
- When an internal standard calibration is used, the IS should be no more than 50% different from the average area of the calibration standards (Section 11.4.2)
- Calibration verification recommended every 10 samples, but only required every 20 samples. Method blanks and control standards are counted in the total. (Section 11.4.1)





Questions?

# Contact Information

Mary Cornell

250 N 17<sup>th</sup> Ave.

Phoenix, AZ 85007

Email: [Mary.Graham@azdhs.gov](mailto:Mary.Graham@azdhs.gov)

Phone: (602) 364-0720