

Arizona Data Qualifiers

Revision 3.0

9/20/2007

Developed by the Sub-committee of the Arizona Environmental Laboratory Advisory Committee

This is an updated list of the Rev. 2.0 Arizona Data Qualifiers dated 11/26/2003, with some new qualifiers added, some obsolete ones deleted and some modified. The new qualifiers are designated in red font. If there was a minor modification to the existing qualifier, it has been highlighted in blue.

Using the following Arizona Data Qualifiers does not automatically denote acceptability to the Regulatory Agency. Arizona Department of Environmental Quality expects that data reported utilizing the following qualifiers, unless stated otherwise, is useable, scientifically valid and defensible. In the laboratory's judgment if the data should not be used for compliance, the T6 qualifier must be used. Other general guidelines for use and application of the following data qualifiers can be found as an attachment to this document (ATTACHMENT A).

Note: For drinking water samples, please use the AZ Drinking Water Data Qualifiers 10/23/2001 [PDF 23K]

Microbiology:

A1 = Too numerous to count.

A2 = Sample incubation period exceeded method requirement.

A3 = Sample incubation period was shorter than method requirement.

A4 = Target organism detected in associated method blank.

A5 = Incubator/water bath temperature was outside method requirements.

A6 = Target organism not detected in associated positive control.

A7 = Micro sample received without adequate headspace.

A8 = Plate count was outside the method's reporting range. Reported value is estimated.

Method/calibration blank:

B1 = Target analyte detected in method blank at or above the method reporting limit.

B2 = Non-target analyte detected in method blank and sample, producing interference.

B3 = Target analyte detected in calibration blank at or above the method reporting limit.

B4 = Target analyte detected in blank at or above method acceptance criteria.

B5 = Target analyte detected in method blank at or above the method reporting limit, but below trigger level or MCL.

B6 = Target analyte detected in calibration blank at or above the method reporting limit, but below trigger level or MCL.

B7 = Target analyte detected in method blank at or above method reporting limit. Concentration found in the sample was 10 times above the concentration found in the method blank.

Confirmation:

C1 = Confirmatory analysis not performed as required by the method.
C2 = deleted
C3 = Qualitative confirmation performed.
C4 = Confirmatory analysis was past holding time.
C5 = Confirmatory analysis was past holding time. Original result not confirmed.
C6 = deleted
C7 = deleted
C8 = Sample RPD between the primary and confirmatory analysis exceeded 40%. Per EPA Method 8000C, the lower value was reported as there was no evidence of chromatographic problems.

Dilution:

D1 = Sample required dilution due to matrix.
D2 = Sample required dilution due to high concentration of target analyte.
D3 = deleted.
D4 = Minimum Reporting Limit (MRL) adjusted to reflect sample amount received and analyzed.
D5 = Minimum Reporting Limit (MRL) adjusted due to sample dilution; analyte was non-detect in the sample.
D6 = Minimum Reporting Limit (MRL) adjusted due to an automatic 10X dilution performed on this sample for the purpose of reporting traditional drinking water analytes for wastewater requirements.

Estimated concentration:

E1 = Concentration estimated. Analyte exceeded calibration range. Reanalysis not possible due to insufficient sample.
E2 = Concentration estimated. Analyte exceeded calibration range. Reanalysis not performed due to sample matrix.
E3 = Concentration estimated. Analyte exceeded calibration range. Reanalysis not performed due to holding time requirements.
E4 = Concentration estimated. Analyte was detected below laboratory minimum reporting limit (MRL).
E5 = Concentration estimated. Analyte was detected below laboratory minimum reporting limit (MRL), but not confirmed by alternate analysis.
E6 = Concentration estimated. Internal standard recoveries did not meet method acceptance criteria.
E7 = Concentration estimated. Internal standard recoveries did not meet laboratory acceptance criteria.

E8 = Analyte reported to MDL per project specification. Target analyte was not detected in the sample.

Hold time:

H1 = Sample analysis performed past holding time.

H2 = Initial analysis within holding time. Reanalysis for the required dilution was past holding time.

H3 = Sample was received and analyzed past holding time.

H4 = Sample was extracted past required extraction holding time, but analyzed within analysis holding time.

H5 = This test is specified to be performed in the field within 15 minutes of sampling; sample was received and analyzed past the regulatory holding time.

BOD/CBOD:

K1 = The sample dilutions set-up for the BOD/CBOD analysis did not meet the oxygen depletion criteria of at least 2 mg/L. Any reported result is an estimated value.

K2 = The sample dilutions set up for the BOD/CBOD analysis did not meet the criteria of a residual dissolved oxygen of at least 1 mg/L. Any reported result is an estimated value.

K3 = deleted.

K4 = deleted.

K5 = The dilution water D.O. depletion was > 0.2 mg/L.

K6 = Glucose/glutamic acid BOD/CBOD was below method acceptance criteria.

K7 = A discrepancy between the BOD and COD results has been verified by reanalysis of the sample for COD.

K8 = Glucose/glutamic acid BOD/CBOD was above method acceptance levels.

Laboratory fortified blank/blank spike:

L1 = The associated blank spike recovery was above laboratory acceptance limits

L2 = The associated blank spike recovery was below laboratory acceptance limits.

L3 = The associated blank spike recovery was above method acceptance limits.

L4 = The associated blank spike recovery was below method acceptance limits.

Matrix spike:

M1 = Matrix spike recovery was high; the associated blank spike recovery was acceptable.

M2 = Matrix spike recovery was low; the associated blank spike recovery was acceptable.

M3 = The spike recovery value is unusable since the analyte concentration in the sample is disproportionate to the spike level. The associated blank spike recovery was acceptable.

M4 = The analysis of the spiked sample required a dilution such that the spike recovery calculation does not provide useful information. The associated blank spike recovery was acceptable.

M5 = Analyte concentration was determined by the method of standard addition (MSA).

M6 = Matrix spike recovery was high. Data reported per ADEQ policy 0154.000.

M7 = Matrix spike recovery was low. Data reported per ADEQ policy 0154.000.

General:

N1 = See case narrative.

N2 = See corrective action report.

N3 = deleted.

N4 = The Minimum Reporting Limit (MRL) verification check did not meet the laboratory acceptance limit.

N5 = The Minimum Reporting Limit (MRL) verification check did not meet the method acceptance limit.

N6 = Data suspect due to quality control failure, reported per data user's request.

Sample quality:

Q1 = Sample integrity was not maintained. See case narrative.

Q2 = Sample received with head space.

Q3 = Sample received with improper chemical preservation.

Q4 = Sample received and analyzed without chemical preservation.

Q5 = Sample received with inadequate chemical preservation, but preserved by the laboratory.

Q6 = Sample was received above recommended temperature.

Q7 = Sample inadequately dechlorinated.

Q8 = Insufficient sample received to meet method QC requirements. Batch QC requirements satisfy ADEQ policies 0154.000 and 0155.000.

Q9 = Insufficient sample received to meet method QC requirements.

Q10 = Sample received in inappropriate sample container.

Q11 = Sample is heterogeneous. Sample homogeneity could not be readily achieved using routine laboratory practices.

Duplicates:

R1 = RPD/RSD exceeded the method acceptance limit.

R2 = RPD/RSD exceeded the laboratory acceptance limit.

R3 = deleted.

R4 = MS/MSD RPD exceeded the method acceptance limit. Recovery met acceptance criteria.

R5 = MS/MSD RPD exceeded the laboratory acceptance limit. Recovery met acceptance criteria.

R6 = LFB/LFBD RPD exceeded the method acceptance limit. Recovery met acceptance criteria.

R7 = LFB/LFBD RPD exceeded the laboratory acceptance limit. Recovery met acceptance criteria.

R8 = Sample RPD exceeded the method acceptance limit.

R9 = Sample RPD exceeded the laboratory acceptance limit.

R10 = deleted.

R11 = The RPD calculation for MS/MSD does not provide useful information due to the varying sample weights when Encore samplers/methanol field preserved samples are used.

Surrogate:

S1 = Surrogate recovery was above laboratory acceptance limits, but within method acceptance limits.

S2 = deleted.

S3 = Surrogate recovery was above laboratory acceptance limits, but within method acceptance limits. No target analytes were detected in the sample.

S4 = Surrogate recovery was above laboratory and method acceptance limits. No target analytes were detected in the sample.

S5 = Surrogate recovery was below laboratory acceptance limits, but within method acceptance limits.

S6 = Surrogate recovery was below laboratory and method acceptance limits.

Reextraction and/or reanalysis confirms low recovery caused by matrix effect.

S7 = Surrogate recovery was below laboratory and method acceptance limits. Unable to confirm matrix effect.

S8 = The analysis of the sample required a dilution such that the surrogate recovery calculation does not provide useful information. The associated blank spike recovery was acceptable.

S9 = deleted.

S10 = Surrogate recovery was above laboratory and method acceptance limits. See case narrative.

S11 = Surrogate recovery was high. Data reported per ADEQ policy 0154.000.

S12 = Surrogate recovery was low. Data reported per ADEQ policy 0154.000.

Method/analyte discrepancies:

T1 = Method approved by EPA, but not yet licensed by ADHS.

T2 = Cited ADHS licensed method does not contain this analyte as part of method compound list.

T3 = Method not promulgated either by EPA or ADHS.

T4 = Tentatively identified compound. Concentration is estimated and based on the closest internal standard.

T5 = Laboratory not licensed for this parameter.

T6 = The reported result cannot be used for compliance purposes.

T7 = Incubator/Oven temperatures were not monitored as required during all days of use.

Calibration verification:

V1 = CCV recovery was above method acceptance limits. This target analyte was not detected in the sample.

V2 = CCV recovery was above method acceptance limits. This target analyte was detected in the sample. The sample could not be reanalyzed due to insufficient sample.

V3 = CCV recovery was above method acceptance limits. This target analyte was detected in the sample, but the sample was not reanalyzed. See case narrative.

V4 = deleted.

V5 = CCV recovery after a group of samples was above acceptance limits. This target analyte was not detected in the sample; acceptable per EPA Method 8000C.

V6 = Data reported from one-point calibration criteria per ADEQ policy 0155.000.

V7 = deleted.

V8 = deleted.

V9 = CCV recovery was below method acceptance limits.

Calibration:

W1= deleted.

W2= deleted.

ATTACHMENT A

“Guidance on the Usage of Data Qualifiers”

These standardized data qualifiers are for use in qualifying analytical results for compliance samples in Arizona to represent events that occurred during analysis. The technical subcommittee has endeavored to develop qualifiers that are succinct and narrow in scope to eliminate broad or multiple interpretations when assessing the impact on data. It must also be noted that due to the specialized nature of the individual qualifiers, it is likely that more than one qualifier may be needed in order to accurately represent the data.

- Note:
1. Using the Arizona Data Qualifiers does not automatically denote acceptability to the Regulatory Agency.
 2. *As specified in the Arizona Adopted Rules, R9-14-615.C.9, for each parameter tested at the laboratory for which quality control acceptance criteria are not specified in the approved method or by EPA or ADEQ,:*
 - a. Use default limits provided in Exhibit II; or*
 - b. Statistically develop limits from historical data*

The laboratory has an option of using ADHS Default Limits which can be accessed at <http://www.azdhs.gov/lab/license/tech/altdefaultlimit.pdf>

Microbiology:

None.

Method/calibration blank:

Apply appropriate qualifier to affected analyte in the blank if target analyte is not detected at \geq RL in the samples. If analytes are detected, then all corresponding analytes for the associated samples should also be qualified.

Confirmation:

For methods that require qualitative confirmation. C3 applies to methods that require quantitative confirmation.

Dilution:

If all analytes are reported from the diluted sample, apply qualifier to the entire sample. Otherwise apply qualifier to each analyte that required dilution.

Estimated concentration:

Appropriate qualifier must be used for any analyte result reported outside the calibration range. Affects data reported outside the calibration range or down to the MDL. E8 is only required if additional clarification is necessary.

Hold time:

Qualify samples appropriately when method extraction and/ or analysis holding time have been exceeded.

BOD/CBOD:

Qualifiers K5, K6, & K8 indicate situations that may impact all results in an analytical run and should be used to qualify all affected samples as well as any affected quality control samples when reported.

Laboratory fortified blank/blank spike:

Appropriate qualifier must be applied to the affected analytes in the Laboratory fortified blank/blank spike and to all corresponding analytes in the associated samples.

Matrix spike:

Appropriate qualifier must be applied to the affected analytes in the matrix spike and should also be added to all corresponding analytes in the associated spiked sample. If a batch spike recovery is outside of the acceptable range, it is permissible to only flag the sample that was spiked and not the other samples in the batch. As required in the Arizona Adopted Rules A.A.C. R9-14-617.8.d, clients must always be informed if the batch QC result is unacceptable whether one of their samples was spiked or not. The laboratory can choose how the unacceptable QC is reported to the client (e.g., cover letter or flag).

The ADEQ policy 0154.000 can be accessed at
<http://www.azdeq.gov/function/business/download/spike8.pdf>

General:

None.

Sample quality:

Flag samples with appropriate qualifier when sample quality may be potentially impacted or when method requirements were not met.

The ADEQ policy 0154.000 can be accessed at

<http://www.azdeq.gov/function/business/download/spike8.pdf>

The ADEQ policy 0155.000 can be accessed at

http://www.azdeq.gov/function/business/download/one_pt3.pdf

Duplicates:

For use with sample, matrix spike, LFB and LFB/blank spike duplicates. Qualify all affected analytes. For MS/MSD or sample duplicates qualify only the original source sample.

Surrogate:

Qualify surrogates appropriately when they do not meet criteria. Surrogate failures in quality control samples will most likely require additional narration. S11 & S12 are used to qualify sample surrogates and only in cases where the Laboratory Fortified Blank/ blank spike has acceptable surrogate recoveries.

Method/analyte discrepancies:

For use with methods or analytes that are not currently approved under the Environmental Laboratory Licensure Rules or for which the lab is not licensed.

Calibration verification:

Appropriate qualifier must be applied to all affected analytes in any samples associated with the calibration verification.

The ADEQ policy 0155.000 can be accessed at

http://www.azdeq.gov/function/business/download/one_pt3.pdf

Calibration:

None.